




EuroMed

Hydrocolloid Health Technologies

Better by Design •

Consumer

Wound Care Products



***Enhanced Healing with
Innovative Hydrocolloid Technology***



Better by Design ●

About EuroMed

EuroMed is a New Jersey based medical device manufacturer with

Expertise in skin adhesive technology

Proprietary and Patented Technology and Process

Professional and Consumer Portfolio

EuroMed products meet rigid safety and certification requirements including ISO 10993 standards and FDA guidelines. EuroMed is also an ISO13485 registered facility as well as FDA QSR compliant.

Our Services

EuroMed is one of a select few custom hydrocolloid manufacturers providing a true turnkey solution.

EuroMed offers contract, private label, and co-marketing services for medical device marketers and distributors worldwide.

EuroMed contract services offer the most extensive, flexible, and unique range of shapes and designs available.

EuroMed offers customized packaging alternatives to meet specific needs.

EuroMed will also store product to provide the convenience of *Just In Time* product delivery to enhance inventory management.

Our contract relationships allow for low cost market entry, complemented by speed-to-market and global market access.

We believe responsiveness to our partners' needs is the foundation for a strong relationship.

Our Quality

We are committed to ensuring an effective quality system that strives to meet customer requirements by aligning our efforts with quality objectives and continuous review of our quality processes.

Hydrocolloid Dressings

Hydrocolloids are substances that convert into a gel when in contact with wound fluid. The moist wound environment created under the dressing speeds the natural healing process, reduces pain and protects against bacterial and other external contaminants.

Hydrocolloid dressings have been clinically shown to

- Facilitate the removal of necrotic tissue and debris through the process of autolytic debridement.

- Reduce the risk of wound infection by providing an effective barrier to external contamination and bacteria.

- Reduces pain and discomfort during dressing removal and application.

- Create an environment that fosters the body's natural healing process resulting in faster healing compared to gauze.

***EuroMed offers a full range of advanced
flexible hydrocolloid dressings in a
variety of shapes and sizes***



Product #: 142510
Size: 21x 42 mm
.82 x 1.6 inches




Product #: 142511
Size: 29 x 59 mm
1.14 x 2.32 inches



Product #: 142512
Size: 17 x 48 mm
.66 x 1.89 inches

***EuroMed offers a full range of advanced
flexible hydrocolloid dressings in a
variety of shapes and sizes***



Product #: 142513
Size: 40 x 64 mm
1.57 x 2.52 inches



Product #: 142514
Size: 23 x 58 mm
.90 x 2.28 inches



Product #: 142515
Size: 30 x 48 mm
1.18 x 1.89 inches

EuroMed also manufactures products for the Professional Market in the following areas

Wound Care

Advanced Hydrocolloid Dressings
Antibacterial Dressings
Dermatological Patches

Ostomy

Ostomy Barriers

Fixation Devices

Male and Female Incontinence
Fixation Devices



EuroMed Inc.
25 Corporate Drive
Orangeburg, NY 10962

P. 1-877-EURO-595
www.euromedinc.com

Better by Design ●