



Better by Design.



Product Catalog
Adhesive Laminates



About EuroMed

EuroMed is a New York-based medical device manufacturer with:

- Expertise in Skin Adhesive Technology
- Proprietary and Patented Technology and Processes
- An Innovative Portfolio of Professional, Consumer, and Roll Good Adhesive Products

EuroMed products meet applicable U.S. and international quality standards. EuroMed is an ISO 13485:2003 and ISO 9001:2000 registered facility, and also is FDA QSR compliant.

Our Services

EuroMed is one of a select few custom hydrocolloid manufacturers providing a true turnkey solution.

EuroMed offers branded (SureSkin® II) and private label products, along with contract and co-marketing services for medical device marketers and distributors worldwide.

EuroMed's contract services offer the most extensive, flexible, and unique range of shapes and designs available.

EuroMed offers customized packaging alternatives to meet specific needs.

EuroMed will store product to provide the convenience of *Just In Time* product delivery to enhance inventory management.

Our contract relationships allow for low-cost market entry, complemented by speed-to-market and global market access.

Our Quality

We are committed to ensuring an effective quality system that strives to meet customer requirements by aligning our efforts with quality objectives and the continuous review of our quality processes.

Hydrocolloid Dressings

Hydrocolloids are substances that convert into a gel when in contact with wound fluid. The moist wound environment created under the dressing speeds the natural healing process, reduces pain, and protects against bacterial and other external contaminants.

Hydrocolloid dressings have been clinically shown to:

- Facilitate the removal of necrotic tissue and debris through the process of autolytic debridement.
- Reduce the risk of wound infection by providing an effective barrier to external contamination and bacteria.
- Reduce pain and discomfort during dressing removal and application.
- Create an environment that fosters the body's natural healing process, resulting in faster healing compared to gauze.





Adhesive Laminates



Non-Sterile Adhesive Laminate Products

EuroMed makes medical grade hydrocolloid adhesive formulations for a variety of professional (medical) and consumer (OTC) needs and indications. The roll stock that EuroMed manufactures follows exacting production protocols and QRS standards in order to provide our customers with the highest level of product quality.

EuroMed's adhesive formulations can be tailored to provide a unique balance between a desired fluid absorption capacity and the strength of adhesion needed for a particular product indication.

For wounds such as pressure ulcers, leg ulcers, or diabetic foot ulcers, EuroMed makes an adhesive formula that is designed to handle increased exudate while providing an adequate level of adhesion to deliver an acceptable wear time.

For wounds such as minor cuts, abrasions, and blisters where wound exudate is minimal and adhesion is of more importance, EuroMed makes an adhesive formula that offers a more suitable alternative that provides increased adhesion capability in relation to a lower absorption capacity.



EuroMed's Adhesive Laminate Formulations

- Adhere to both dry and moist skin surfaces.
- Deliver gentle and secure adhesion that does not irritate the skin.
- Have fluid absorption capacity that is adjustable based on adhesive formula.
- Provide barrier properties to external contaminants with secure adhesion to skin.
- Are flexible and conformable and mold to difficult anatomical contours.
- Come with silicone release liners to facilitate removal from adhesive and application to skin.
- Contain no animal derived byproducts.
- Are latex free.

Roll good laminate product Data Sheets are available upon request.

Common Indications for Use

- Professional wound dressings for pressure, leg, and diabetic ulcers.
- Post-op dressings for suture line protection.
- Peripheral I.V. and central line management.
- Ostomy adhesive barriers to protect the skin from stomal effluent.
- Incontinence fixation adhesive to provide secure adhesion of external catheters.
- Consumer first aid dressings for minor cuts, abrasions, and burns.
- Consumer blemish dressings.
- Consumer foot and heel dressings.



Better by Design.



The information presented here is intended as background data to assist in product evaluation and verification. The product features, properties, characteristics, and performance data have been determined through sampling and testing during product manufacture. A material used in a particular application must be tested by the customer for suitability for the intended use, and the customer assumes all risk and liability whatsoever in connection therewith. This document will not form any part of ensuing sales contracts.

U.S. Patents: 6309500, 6143798, 5998694, 5571080, 5704905
Other International patents apply.
Additional patents pending.

Sureskin is a registered trademark of Euromed, Inc. in Austria, Sweden, Spain, Denmark, Germany and France.
Euromed, Inc. has pending trademark applications in other foreign jurisdictions.

EuroMed Inc.
25 Corporate Drive, Orangeburg, NY 10962

info@euromedinc.com | P: 845-359-4039 • 877-Euro-595 | F: 845-359-1315
www.euromedinc.com